

NO: CE3002

## **SYBR Green/SYBR Premix Ex Taq Kit**

### **Introduction**

SYBR Premix Ex Taq Kit is used for intercalator based real time PCR using SYBR®Green. Rapid reaction mixtures are based on a 2×concentration premixed with SYBR®Green I at a concentration appropriate for real time monitoring.

A combination of Ex Taq™ HS, a hot start PCR enzyme that uses an anti-Taq antibody, and a buffer optimized for real time PCR allows high amplification efficiency and high detection sensitivity in real time PCR. The 2×premixed reagent also contains Tli RNaseH, a heat-resistant RNase H, which minimizes inhibition of PCR due to residual mRNA when using cDNA as template. This product is suitable for high-speed PCR and allows accurate assay and detection of targets, making it possible to conduct real time PCR analyses with good reproducibility and high reliability.

### **Storage**

2-8°C avoid light for 6 months, -20°C for two years. Do not freeze (Return to room temperature before use).

### **Components**

1. SYBR Premix Ex Taq (2×Conc.) (5units,1.0 mL)
2. Specification: 500preps.

### **DNA Preparation--Please select ECALBIO of DNA Purification kit (Spin Column)**

#### ***-Blood***

1. Take 500μL anticoagulant blood to a 1.5 mL clean centrifuge tube, add 1 ml of ddH<sub>2</sub>O, shake for 30s .
2. Centrifuge at 8,000 rpm for 5 min, discard the supernatant.
3. Repeat washing until no red precipitate.
4. Add 1mL physiological saline, shake 15s, centrifuge at 10,000 rpm for 5 min, discard the supernatant.
5. Add 50μL DNA extraction and mix with precipitation.
6. 100°C 10 min. Centrifuged at 13,000 rpm for 3min.
7. Take the supernatant 4μL for PCR reactions.

#### ***-Animal tissue, food or feed***

1. Take about 30mg homogenized samples, and put into a 1.5mL clean centrifuge tube, add 100μL DNA extraction and mix well (the DNA extraction must dissolve at room temperature and fully mixed before use).
2. Water bath at 56°C for 30 min, then at 100°C for 10 min.
3. Centrifuge at 13,000rpm for 5 min.
4. Transfer the supernatant 5μL for PCR reactions.

### **Real Time PCR Reaction**

1. Take out qPCR MIX from the kit, melt it at room temperature and oscillate it, then mix it with 10000 rpm centrifugal 10s;
2. Each of the 25μL system, the upstream and downstream primers were 0.5μL, add 12.5μL SYBR Green Premix and 9μL ddH<sub>2</sub>O respectively, is up to 25μL system, mixed and amplified.

PCR Solution	PCR System (20μL)	PCR System (25μL)
DNA Sample	2μL	2.5μL
Upstream primers	0.4μL	0.5μL
Downstream primers	0.4μL	0.5μL
2×SYBR Green Premix	10μL	12.5μL
ddH <sub>2</sub> O	7.20μL	9.0μL

3. Example: IQTM5, ABI PRISM7500 and other instruments using thin walled tubes: circulation conditions: 95°C for 3 minutes, then 95°C for 10-15s, 60°C for 30s, 72°C for 30s, 40 cycles. The user should set the corresponding amplification conditions according to their own primers, and the above conditions are for reference only.

PCR amplification	40 Cycle	Step 1	95 °C for 3 min
		step 2	94 °C for 10-15 s
			60 °C for 20 s
			72 °C for 30 s

### Results determination

- Under the test conditions established, Ct value  $\leq 10$ , dilute the DNA Sample, the real time PCR test need to be done again.
- 20 < Ct value  $\leq 30$ , quantitative analysis is accuracy.
- If 30 < Ct value  $\leq 35$ , need to increase the volume of DNA sample, improve PCR amplification to get correct result.
- Ct value > 35, check amplification curve. If the amplification curve is logarithmic amplification curve, it was suspected positive, otherwise judged as negative.